

Entermed patient chair

Fully automatic patient examination chair



The fully automatic examination chair is comfortable for both doctor and patient. It can be rotated 350 degrees and has controls and user definable memory buttons on the side of the backrest for easy accessibility. It comes with foot control and can be controlled by the Futurent III ENT treatment unit. The memory buttons can be programmed quickly and simply. The armrests are hinged and the footrest can be folded out of the way when not required. Its headrest is double articulated and height adjustable to allow the doctor to place the patient's head in the desired position. The chair is very stable in all positions due to its clever design. It can be supplied with retractable wheels.

*The memory buttons can be programmed
quickly and simply*

Benefits at a glance

- Elegant design
- Synchronised backrest and legrest
- Stainless steel base
- User definable memory buttons
- Rotation lock on both sides
- Large variety of colours and upholstery design
- Hinged armrests
- Folding footrest
- Electric height and backrest adjustment
- Dual-function foot control (up-down, front-back)
- Foot control and set of controls on backrest
- Available with retractable wheels



Horizontal position up to -20°



Double-articulated headrest

Technical specifications

Base dimensions (L x D)	56 x 71 cm
Seat dimensions (L x D)	51 x 54 cm
Maximum seat height	79 cm (without headrest)
Minimum seat height	54 cm (without headrest)
Overall length	147 cm (without headrest & footrest)
Rotation	350°
Maximum load	150 kg
Vertical position	up to +5°
Horizontal position	up to -20°

Product code: 760110

Entermed patient chair
Fully automatic patient examination chair

Product code: 760110-WH

Retractable wheels